

Dual Layer Amniotic Membrane

- ❑ Stores dry at room temperature
- ❑ Suitable for immediate use off the shelf
- ❑ Terminally sterilized with gamma radiation to SAL 10⁻⁶ in accordance with ISO 11137
- ❑ Not chemically cross-linked
- ❑ Three year shelf life
- ❑ Axograft™ Amniotic Membrane (Dual Layer) version processed with surface layer dimpling for ease of graft placement / handling
- ❑ Collection of the donor tissue is performed in an aseptic manner by appropriately licensed tissue establishments
- ❑ Placental donors go through a rigorous pre-screening qualification
- ❑ Processed in accordance with the safety guidelines provided by the U.S. Food and Drug Administration (FDA) - Human Cellular and Tissue-based Products (HCT/P) (21 CFR Part 1271)
- ❑ Processed in accordance with the standards from the American Association of Tissue Banks (AATB)



Part #	Description
X090-0025-AMN22CMCC	2 x 2cm Amniotic Membrane, Axograft
X090-0025-AMN44CMCC	4 x 4cm Amniotic Membrane, Axograft
X090-0025-AMN46CMCC	4 x 6cm Amniotic Membrane, Axograft
X090-0025-AMN48CMCC	4 x 8cm Amniotic Membrane, Axograft

INDICATIONS: See Package Insert for a more complete listing of indications, contraindications, warnings, precautions, and other important information.

LIMITED WARRANTY and DISCLAIMER: X-spine products have a limited warranty against defects and workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are disclaimed.

WARNING: In the USA, this product has labeling limitations. See package insert for complete information.

CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

X-spine®, the X-spine logo, the Axograft™ logo, and Axograft™ mark are all trademarks or service marks of X-spine Systems, Inc.

The "Axograft" name is a registered trademark of X-spine Systems, Inc.; Amniotic Membrane supplied by Bone Bank Allografts.

Patents Pending.

All products are not currently available in all markets. © 2016 X-spine Systems, Inc., All rights reserved.

