


# Precision-Milled Allografts



BacFast HD



BacFast HD Surface via a Scanning Electron Microscope (SEM) at a magnification of 22x

## BacFast HD

BacFast HD is a facet stabilization dowel designed with a focus on osteoconductivity and preventing graft migration, making it an ideal graft to be used as a complement to common spinal fusion procedures as it preserves the posterior spine. Utilizing Bacterin's proprietary demineralization technology, BacFast HD is hyper-demineralized to expose the growth factors and BMPs inherent to cortical bone. With the benefits of HD technology and increased collagen surface area, BacFast HD also provides the graft with osteoinductive properties without compromising the structural integrity of the graft.

## Precision-Milled Allografts



BacFast HD on Inserter  
Instrumentation Surgical Kits  
available upon request

- |         |                                  |
|---------|----------------------------------|
| 100511  | BacFast HD Facet Dowel, 5 mm     |
| 100512  | BacFast HD Facet Dowel, 6 mm     |
| 1563305 | Cervical Spacer, Lordotic, 5 mm  |
| 1563306 | Cervical Spacer, Lordotic, 6 mm  |
| 1563307 | Cervical Spacer, Lordotic, 7 mm  |
| 1563308 | Cervical Spacer, Lordotic, 8 mm  |
| 1563309 | Cervical Spacer, Lordotic, 9 mm  |
| 1563310 | Cervical Spacer, Lordotic, 10 mm |