

CALIX[™]
(CERVICAL INTERBODY SPACER)

CALIX·A[™]
(PEEK LUMBAR SYSTEM)




CALIX·T[™]
(PEEK LUMBAR SYSTEM)

CALIX·P[™]
(PEEK LUMBAR SYSTEM)



CALIX™

CERVICAL INTERBODY SPACER



Convex



Lordotic







● 5 mm ● 6 mm ● 7 mm ● 8 mm ● 9 mm ● 10 mm ● 11 mm ● 12 mm



CALIX·A™

PEEK LUMBAR SYSTEM









☒ Large grafting area

☒ Aggressive teeth help prevent migration

☒ Bullet Nose for ease of insertion

☒ Tantalum markers for implant visualization





Disk Preparation Instruments




CALIX™

CERVICAL INTERBODY SPACER


Implants

Item #	Description
X003-0005	11mm x 14mm x 5mm x 7°, Lordotic
X003-0006	11mm x 14mm x 6mm x 7°, Lordotic
X003-0007	11mm x 14mm x 7mm x 7°, Lordotic
X003-0008	11mm x 14mm x 8mm x 7°, Lordotic
X003-0009	11mm x 14mm x 9mm x 7°, Lordotic
X003-0010	11mm x 14mm x 10mm x 7°, Lordotic
X003-0011	11mm x 14mm x 11mm x 7°, Lordotic
X003-0012	11mm x 14mm x 12mm x 7°, Lordotic
X003-0105	11mm x 14mm x 5mm Convex
X003-0106	11mm x 14mm x 6mm Convex
X003-0107	11mm x 14mm x 7mm Convex
X003-0108	11mm x 14mm x 8mm Convex
X003-0109	11mm x 14mm x 9mm Convex
X003-0110	11mm x 14mm x 10mm Convex
X003-0111	11mm x 14mm x 11mm Convex
X003-0112	11mm x 14mm x 12mm Convex



CALIX·A™

(PEEK LUMBAR SYSTEM)



Item #	Description
X034-0320	36mm x 24mm x 8mm x 8°
X034-0321	36mm x 24mm x 10mm x 8°
X034-0322	36mm x 24mm x 12mm x 8°
X034-0323	36mm x 24mm x 14mm x 8°
X034-0324	36mm x 24mm x 16mm x 8°
X034-0325	36mm x 24mm x 18mm x 8°
X034-0326	36mm x 24mm x 20mm x 8°
X034-0330	36mm x 24mm x 8mm x 12°
X034-0331	36mm x 24mm x 10mm x 12°
X034-0332	36mm x 24mm x 12mm x 12°
X034-0333	36mm x 24mm x 14mm x 12°
X034-0334	36mm x 24mm x 16mm x 12°
X034-0335	36mm x 24mm x 18mm x 12°
X034-0336	36mm x 24mm x 20mm x 12°

CALIX·P™

(PEEK LUMBAR SYSTEM)



Item #	Description
X034-0171	22mm x 10mm x 6mm
X034-0172	22mm x 10mm x 7mm
X034-0173	22mm x 10mm x 8mm
X034-0174	22mm x 10mm x 9mm
X034-0175	22mm x 10mm x 10mm
X034-0176	22mm x 10mm x 11mm
X034-0177	22mm x 10mm x 12mm
X034-0178	22mm x 10mm x 13mm
X034-0179	22mm x 10mm x 14mm
X034-0180	22mm x 10mm x 15mm
X034-0181	22mm x 10mm x 16mm
X034-0240	26mm x 10mm x 6mm
X034-0241	26mm x 10mm x 8mm
X034-0246	26mm x 10mm x 9mm
X034-0242	26mm x 10mm x 10mm
X034-0247	26mm x 10mm x 11mm
X034-0243	26mm x 10mm x 12mm
X034-0248	26mm x 10mm x 13mm
X034-0244	26mm x 10mm x 14mm
X034-0387	26mm x 10mm x 6mm x 5°
X034-0389	26mm x 10mm x 8mm x 5°
X034-0390	26mm x 10mm x 9mm x 5°
X034-0391	26mm x 10mm x 10mm x 5°
X034-0392	26mm x 10mm x 11mm x 5°
X034-0393	26mm x 10mm x 12mm x 5°
X034-0394	26mm x 10mm x 13mm x 5°
X034-0395	26mm x 10mm x 14mm x 5°
X034-0396	26mm x 10mm x 15mm x 5°

CALIX·T™

(PEEK LUMBAR SYSTEM)



Item #	Description
X034-0280	28mm x 10mm x 6mm
X034-0281	28mm x 10mm x 7mm
X034-0282	28mm x 10mm x 8mm
X034-0283	28mm x 10mm x 9mm
X034-0284	28mm x 10mm x 10mm
X034-0285	28mm x 10mm x 11mm
X034-0286	28mm x 10mm x 12mm
X034-0287	28mm x 10mm x 13mm
X034-0288	28mm x 10mm x 14mm
X034-0289	28mm x 10mm x 15mm
X034-0290	28mm x 10mm x 16mm

WARNING: In the USA, this product has labeling limitations. See package insert for complete information.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

X-spineSM the X-spine logoSM and CalixTM are trademarks of X-spine Systems, Inc.

Products Patented and Patents Pending

All products are not currently available in all markets.

© 2012 X-spine Systems, Inc., All rights reserved.

 X-spineSM
X-treme Innovations

X-spine Systems, Inc.

452 Alexandersville Rd., Miamisburg, OH 45342

Phone: 800-903-0640 • Direct: 937-847-8400 • Fax: 937-847-8410

www.x-spine.com



F-1000.17 Rev A I2M/IC
EC 60-1253

