


# FORTEX<sup>®</sup>

Pedicle Screw System

Simplicity  
Reinvented



- ❑ Acute thread angle reduces the risk of backout and migration of the screw after implantation
- ❑ Rods available in both Titanium and Cobalt Chromium
- ❑ 70° of polyaxial motion

 X·spine<sup>®</sup>


An Xtant Medical Company



## Fortex Pedicle Screws and Cap



Item#	Description
X022-0080	Fortex, 4.75mm x 30mm Cannulated Screw
X022-0081	Fortex, 4.75mm x 35mm Cannulated Screw
X022-0082	Fortex, 4.75mm x 40mm Cannulated Screw
X022-0083	Fortex, 4.75mm x 45mm Cannulated Screw
X022-0084	Fortex, 4.75mm x 50mm Cannulated Screw
X022-0085	Fortex, 4.75mm x 55mm Cannulated Screw
X022-0086	Fortex, 5.5mm x 30mm Cannulated Screw
X022-0087	Fortex, 5.5mm x 35mm Cannulated Screw
X022-0088	Fortex, 5.5mm x 40mm Cannulated Screw
X022-0089	Fortex, 5.5mm x 45mm Cannulated Screw
X022-0090	Fortex, 5.5mm x 50mm Cannulated Screw
X022-0091	Fortex, 5.5mm x 55mm Cannulated Screw
X022-0092	Fortex, 6.5mm x 30mm Cannulated Screw
X022-0093	Fortex, 6.5mm x 35mm Cannulated Screw
X022-0094	Fortex, 6.5mm x 40mm Cannulated Screw
X022-0095	Fortex, 6.5mm x 45mm Cannulated Screw
X022-0096	Fortex, 6.5mm x 50mm Cannulated Screw
X022-0097	Fortex, 6.5mm x 55mm Cannulated Screw
X022-0098	Fortex, 7.5mm x 30mm Cannulated Screw
X022-0099	Fortex, 7.5mm x 35mm Cannulated Screw
X022-0100	Fortex, 7.5mm x 40mm Cannulated Screw
X022-0101	Fortex, 7.5mm x 45mm Cannulated Screw
X022-0102	Fortex, 7.5mm x 50mm Cannulated Screw
X022-0103	Fortex, 7.5mm x 55mm Cannulated Screw
X022-0104	Fortex, 8.25mm x 30mm Cannulated Screw
X022-0105	Fortex, 8.25mm x 35mm Cannulated Screw
X022-0106	Fortex, 8.25mm x 40mm Cannulated Screw
X022-0107	Fortex, 8.25mm x 45mm Cannulated Screw
X022-0108	Fortex, 8.25mm x 50mm Cannulated Screw
X022-0109	Fortex, 8.25mm x 55mm Cannulated Screw
X022-0042	Fortex Cap



## Fortex Cross Connectors

Item#	Description
X070-0010	35 - 45mm Cross Bar Connector
X070-0011	45 - 60mm Cross Bar Connector
X070-0012	60 - 80mm Cross Bar Connector

## Fortex Self Drilling Pedicle Screws\* and Cap



Item#	Description
X022-5530SD	Fortex, 5.5mm x 30mm Self-Drilling Screw
X022-5535SD	Fortex, 5.5mm x 35mm Self-Drilling Screw
X022-5540SD	Fortex, 5.5mm x 40mm Self-Drilling Screw
X022-5545SD	Fortex, 5.5mm x 45mm Self-Drilling Screw
X022-5550SD	Fortex, 5.5mm x 50mm Self-Drilling Screw
X022-5555SD	Fortex, 5.5mm x 55mm Self-Drilling Screw
X022-6530SD	Fortex, 6.5mm x 30mm Self-Drilling Screw
X022-6535SD	Fortex, 6.5mm x 35mm Self-Drilling Screw
X022-6540SD	Fortex, 6.5mm x 40mm Self-Drilling Screw
X022-6545SD	Fortex, 6.5mm x 45mm Self-Drilling Screw
X022-6550SD	Fortex, 6.5mm x 50mm Self-Drilling Screw
X022-6555SD	Fortex, 6.5mm x 55mm Self-Drilling Screw
X022-7530SD	Fortex, 7.5mm x 30mm Self-Drilling Screw
X022-7535SD	Fortex, 7.5mm x 35mm Self-Drilling Screw
X022-7540SD	Fortex, 7.5mm x 40mm Self-Drilling Screw
X022-7545SD	Fortex, 7.5mm x 45mm Self-Drilling Screw
X022-7550SD	Fortex, 7.5mm x 50mm Self-Drilling Screw
X022-7555SD	Fortex, 7.5mm x 55mm Self-Drilling Screw
X022-0042	Fortex Cap

\*Special Order

## Fortex Rods

Item#	Description
N60001073	35mm Straight Rod
N60000300	40mm Straight Rod
N60001074	45mm Straight Rod
N60001075	50mm Straight Rod
N60001076	55mm Straight Rod
N60000301	60mm Straight Rod
N60001077	65mm Straight Rod
N60001060	70mm Straight Rod
N60000302	80mm Straight Rod
N60001080	90mm Straight Rod
N60000303	100mm Straight Rod
N60000304	120mm Straight Rod
N60000307	180mm Straight Rod
N60000309	300mm Straight Rod
X077-0420	5.5mm Cobalt Chromium Straight Rod, 500mm*
N60001091	35mm Continuous Radius Rod
N60001092	40mm Continuous Radius Rod
N60001093	45mm Continuous Radius Rod
N60001068	50mm Continuous Radius Rod
N60001069	60mm Continuous Radius Rod
N60001070	70mm Continuous Radius Rod
N60001071	80mm Continuous Radius Rod
N60001098	90mm Continuous Radius Rod
N60001072	100mm Continuous Radius Rod

\*Special Order



## X-spine Systems, Inc.

452 Alexandersville Rd., Miamisburg, OH 45342  
 Phone: 800-903-0640 • Direct: 937-847-8400 • Fax: 937-847-8410  
[www.x-spine.com](http://www.x-spine.com)

INDICATIONS: See Package Insert for a more complete listing of indications, contraindications, warnings, precautions, and other important information.

LIMITED WARRANTY and DISCLAIMER: X-spine products have a limited warranty against defects and workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are disclaimed.

WARNING: In the USA, this product has labeling limitations. See package insert for complete information. CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

X-spine®, the X-spine logo®, and Fortex® are trademarks or servicemarks of X-spine Systems, Inc.

Patents Pending.

© 2016 X-spine Systems, Inc., All rights reserved.