



IRIX-A™


Lumbar Integrated Fusion System





- ❑ Integrated titanium ring, surrounded by outer PEEK ring, promotes increased torsional strength and shear performance
- ❑ Resilient locking mechanism, built upon proven lumbar fusion technology
- ❑ Frictional Titanium Plasma Coating on bone contact areas
- ❑ Implants are sterile packaged for product quality




 Titanium Plasma Coated (PC)
 PEEK provides fluoroscopic
 endplate visualization


 Locking arm provides
 tactile, visual and
 audible feedback




 Screw orientation 35°
 cephalad/caudal with
 up to 6° conical screw
 angulation

 Titanium ring provides
 better x-ray visibility


Irix-A™
33mm x 26mm Implants



Irix-A™
38mm x 28mm Implants



Irix-A™
42mm x 30mm Implants



Item #	Description
X080-332611-08PC-STR	33 X 26 X 8°, 11mm, PC (Sterile)
X080-332613-08PC-STR	33 X 26 X 8°, 13mm, PC (Sterile)
X080-332615-08PC-STR	33 X 26 X 8°, 15mm, PC (Sterile)
X080-332617-08PC-STR	33 X 26 X 8°, 17mm, PC (Sterile)
X080-332619-08PC-STR	33 X 26 X 8°, 19mm, PC (Sterile)
X080-332611-12PC-STR	33 X 26 X 12°, 11mm, PC (Sterile)
X080-332613-12PC-STR	33 X 26 X 12°, 13mm, PC (Sterile)
X080-332615-12PC-STR	33 X 26 X 12°, 15mm, PC (Sterile)
X080-332617-12PC-STR	33 X 26 X 12°, 17mm, PC (Sterile)
X080-332619-12PC-STR	33 X 26 X 12°, 19mm, PC (Sterile)

Item #	Description
X080-382811-08PC-STR	38 X 28 X 8°, 11mm, PC (Sterile)
X080-382813-08PC-STR	38 X 28 X 8°, 13mm, PC (Sterile)
X080-382815-08PC-STR	38 X 28 X 8°, 15mm, PC (Sterile)
X080-382817-08PC-STR	38 X 28 X 8°, 17mm, PC (Sterile)
X080-382819-08PC-STR	38 X 28 X 8°, 19mm, PC (Sterile)
X080-382811-12PC-STR	38 X 28 X 12°, 11mm, PC (Sterile)
X080-382813-12PC-STR	38 X 28 X 12°, 13mm, PC (Sterile)
X080-382815-12PC-STR	38 X 28 X 12°, 15mm, PC (Sterile)
X080-382817-12PC-STR	38 X 28 X 12°, 17mm, PC (Sterile)
X080-382819-12PC-STR	38 X 28 X 12°, 19mm, PC (Sterile)

Item #	Description
X080-423011-08PC-STR	42 X 30 X 8°, 11mm, PC (Sterile)
X080-423013-08PC-STR	42 X 30 X 8°, 13mm, PC (Sterile)
X080-423015-08PC-STR	42 X 30 X 8°, 15mm, PC (Sterile)
X080-423017-08PC-STR	42 X 30 X 8°, 17mm, PC (Sterile)
X080-423019-08PC-STR	42 X 30 X 8°, 19mm, PC (Sterile)
X080-423011-12PC-STR	42 X 30 X 12°, 11mm, PC (Sterile)
X080-423013-12PC-STR	42 X 30 X 12°, 13mm, PC (Sterile)
X080-423015-12PC-STR	42 X 30 X 12°, 15mm, PC (Sterile)
X080-423017-12PC-STR	42 X 30 X 12°, 17mm, PC (Sterile)
X080-423019-12PC-STR	42 X 30 X 12°, 19mm, PC (Sterile)

 **Self-Drilling Screws**


Item #	Description
X080-5525SD	5.5 X 25mm, Self-Drilling Screw
X080-5530SD	5.5 X 30mm, Self-Drilling Screw
X080-5535SD	5.5 X 35mm, Self-Drilling Screw
X080-5540SD	5.5 X 40mm, Self-Drilling Screw

 **Self-Tapping Screws**


Item #	Description
X080-5525ST	5.5 X 25mm, Self-Tapping Screw
X080-5530ST	5.5 X 30mm, Self-Tapping Screw
X080-5535ST	5.5 X 35mm, Self-Tapping Screw
X080-5540ST	5.5 X 40mm, Self-Tapping Screw

 **Self-Tapping Rescue Screws**


Item #	Description
X080-6025	6.0 X 25mm, Self-Tapping Rescue Screw
X080-6030	6.0 X 30mm, Self-Tapping Rescue Screw
X080-6035	6.0 X 35mm, Self-Tapping Rescue Screw
X080-6040	6.0 X 40mm, Self-Tapping Rescue Screw



X-spine Systems, Inc.
 452 Alexandersville Rd., Miamisburg, OH 45342
 Phone: 800-903-0640 • Direct: 937-847-8400 • Fax: 937-847-8410
www.x-spine.com

INDICATIONS: See Package Insert for a more complete listing of indications, contraindications, warnings, precautions, and other important information.

LIMITED WARRANTY and DISCLAIMER: X-spine products have a limited warranty against defects and workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are disclaimed.

WARNING: In the USA, this product has labeling limitations. See package insert for complete information. CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

X-spineSM, the X-spine logo, and Irix-ATM are trademarks or servicemarks of X-spine Systems, Inc.

Patents Pending.

© 2016 X-spine Systems, Inc., All rights reserved.

F-1000.43 Rev B 2/16

