


Spider[®]

Cervical Plating System

Revolutionary Locking: Audible, Tactile & Visual



- ❑ Resilient locking arm technology
- ❑ Serrated bottom surface for greater grip
- ❑ Monoblock construction




An Xtant Medical Company



Locking arms move aside to let screw pass and then “click” back over to capture and lock screw.

Resilient locking arms provide single step locking with 3 forms of locking feedback: audible, tactile, and visual.




Spider® Plates

Item Number	Description
N60000100	1 Level 22mm Plate
N60000101	1 Level 24mm Plate
N60000102	1 Level 26mm Plate
N60000103	1 Level 28mm Plate
N60000104	1 Level 30mm Plate
N60000105	1 Level 32mm Plate
N60000106	1 Level 34mm Plate
N60000107	2 Level 36mm Plate
N60000108	2 Level 39mm Plate
N60000109	2 Level 42mm Plate
N60000110	2 Level 45mm Plate
N60000111	2 Level 48mm Plate
N60000112	2 Level 51mm Plate
N60000113	2 Level 54mm Plate
N60000114	3 Level 53mm Plate
N60000115	3 Level 56mm Plate
N60000116	3 Level 59mm Plate
N60000117	3 Level 62mm Plate
N60000118	3 Level 65mm Plate
N60000129	3 Level 68mm Plate
N60000119	3 Level 71mm Plate
N60000120	3 Level 74mm Plate
N60000121	3 Level 77mm Plate
N60000122	4 Level 68mm Plate
N60000123	4 Level 72mm Plate
N60000124	4 Level 76mm Plate
N60000125	4 Level 80mm Plate
N60000126	4 Level 84mm Plate
N60000127	4 Level 88mm Plate
N60000128	4 Level 92mm Plate



Spider® Screws

Item Number	Description
N60000135	4mm x 12mm Variable Screw
N60000136	4mm x 14mm Variable Screw
N60000137	4mm x 16mm Variable Screw
N60000140	4.25mm x 12mm Variable Screw
N60000141	4.25mm x 14mm Variable Screw
N60000142	4.25mm x 16mm Variable Screw
N60000145	4mm x 12mm Fixed Screw
N60000146	4mm x 14mm Fixed Screw
N60000147	4mm x 16mm Fixed Screw
N60000150	4.25mm x 12mm Fixed Screw
N60000151	4.25mm x 14mm Fixed Screw
N60000152	4.25mm x 16mm Fixed Screw



*Four level plates available as special order.



X-spine Systems, Inc.

452 Alexandersville Rd., Miamisburg, OH 45342
 Phone: 800-903-0640 • Direct: 937-847-8400 • Fax: 937-847-8410
www.x-spine.com

INDICATIONS: See Package Insert for a more complete listing of indications, contraindications, warnings, precautions, and other important information.

LIMITED WARRANTY and DISCLAIMER: X-spine products have a limited warranty against defects and workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are disclaimed.

WARNING: In the USA, this product has labeling limitations. See package insert for complete information. CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

X-spine®, the X-spine logo®, and Spider® are trademarks or servicemarks of X-spine Systems, Inc.

Patents Pending.

© 2016 X-spine Systems, Inc., All rights reserved.

F-1000.9 Rev C 6/16

